

Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that our product

smoov (one) O10

including the lithium-ion batteries meets the requirements of

UN Manual of Tests and Criteria" Part III subsection 38.3

and therefore are suitable for transportation by passengers according

IATA DGR Chapter 2.3.2.4

The rated energy of the built-in drive unit battery (which is not intended for removal by the user) is **216 Wh**. The rated energy of the control unit (item number 1594326) battery (also not removable by the user) is 9,9 Wh.

Since drive unit and control unit are portable and contain batteries, we recommend that you take them with you to the passenger cabin. Both units must be switched off before boarding and protected from damage and unintentional activation. The pilot in command must be informed about the location of the mobility aid.

We strongly recommend contacting the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defective devices or batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2024**.

Albstadt, 13th of November 2023

Alber GmbH



Johannes Beltman
Managing Director