

Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that our product

e-motion M25 (incl. DuoDrive)

including the lithium-ion batteries fulfil the requirements of

IATA chapter 2.3.2.4

and are therefore suitable for transportation by passengers.

The pilot in command must be informed about the location of the mobility aid.

The rated energy of the built-in wheel drive unit battery (not intended for removal by the user) is < **160 Wh**. The rated energy of the DuoDrive control unit is < **11 Wh**.

The e-motion M25 can be set into a dedicated "flight mode". In this mode the system is completely deactivated. While in this mode the e-motion M25 can only be switched on by a coded procedure. Unintended switching on by luggage is not possible. Furthermore, the Bluetooth device of the system is completely switched off. As the DuoDrive control unit is portable and contains a battery, we recommend that you take it with you into the passenger cabin.

We strongly recommend contacting the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2024**.

Albstadt, 04th of November 2024

Alber GmbH



Joost Beltman
Managing Director